



Pavia Rectal Temp Short Flex Tip

Fast · Accurate · Comfortable

- Accurate results in 6 seconds
- Short Flexible Tip (probe) provides comfort for the pet
- Length of probe just 2 inches including metal end
- Handheld cordless operation
- Operates on 2 AAA batteries
- Weighs just 1 oz, with batteries
- Four digit LCD display
- Low battery indicator
- Last temperature is recorded in memory
- Manual or automatic 60 second shut-off
- Temperature range 89.6 degrees F-107.6 degrees F
- 90 day warranty

For use on dogs, cats, and other small animals only.

Average life- approximately 9 months depending on how many temps you take and clinic breakage.

PAVIA PROBE COVERS

The Pavia Rectal Temp Short Flex Tip requires the use of probe covers for sanitary reasons, due to the flexible, porous, rubber material of the probe. Pavia Probe Covers are designed to fit the Pavia Short Flex Tip thermometer. Sold Separately. 100 covers per box.



Pavia Sales Group, Inc.
3505 Ranier Lane North • Plymouth, MN 55447
Phone: 866-593-0790 • Fax: 763-473-1835
www.paviasalesgroup.com

